


Protocol & Procedure for Placement of the OCO Biomedical Engage™ Two-Piece Implant



Proper Drill Sequence for Engage™ Two-Stage Implant (3.25, 4.0, 5.0, 6.0, 7.0 & 8.0mm diameters)

- 1 - #8 High-speed Surgical Bur** - Penetrate soft tissue if flapless, otherwise use to divot crest
- 2 - 1.8mm Pilot Drill** - Place drill stop onto pilot drill and drill full length stopping at soft tissue
- 3 - Paralleling Pin** - Place paralleling pin to evaluate implant size and trajectory
- 4 - Tissue Punch** - Remove soft tissue if placing implant flapless
- 5 - 1.8mm Pilot Drill** - Take pilot drill to full length of intended implant
- 6 - Final Drill** - Take final drill to full length of osteotomy. Use upper drill for maxillary/soft bone and lower for mandibular/hard bone.*
- 8 - Internal Hex Driver** - Drive implant into final position

*6.0mm, 7.0mm and 8.0mm diameter Engage™ implants must have osteotomy created by stepping up in serial progression

For questions on the Engage™ implant placement and restorative techniques please visit www.ocobiomedical.com, call 800-228-0477, +1 (505)-293-0025 international or email sales@ocobiomedical.com.