


Protocol & Procedure for Placement of the OCO Biomedical SDI 2.9mm and 3.0 One-Piece Implant



Proper Drill Sequence for 2.9 SDI & 3.0 One-Piece Implant

- 1 - #8 High-speed Surgical Bur - Penetrate soft tissue if flapless, otherwise use to divot crest
- 2 - 1.8mm Pilot Drill - Place 8mm drill stop onto pilot drill and drill full length stopping at soft tissue
- 3 - Paralleling Pin - Place paralleling pin to evaluate implant size and trajectory
- 4 - Tissue Punch - Remove soft tissue if placing implant flapless
- 5 - Countersink Drill - Take top of blades approximately 1 - 1.5mm below soft tissue
- 6 - 1.8mm Pilot Drill - Take pilot drill to full length of intended implant
- 7 - Final Drill - Take final drill to full length of osteotomy
- 8 - Driver: C & B or O-Ball Driver - Drive implant into final position

For questions on the 2.9(mm) SDI or 3.0(mm) implant placement and restorative techniques please visit www.ocobiomedical.com, call 800-228-0477, +1 (505)-293-0025 international or email sales@ocobiomedical.com.

OCO Biomedical
9550 San Mateo Blvd. NE, Suite C
Albuquerque, NM 87113
www.ocobiomedical.com

Toll Free: (800) 228-0477
International: + 1(505) 293-0025
Direct: (505) 293-0025
Fax: (505) 293-0447

Designed and Manufactured in USA
U.S. Patent: US 8,277,218,B2
ISO 13485:2003 CERTIFIED

